

Quick Reference Guide – AUSCARE Episodes View

Investigations **Ordering**

Sign Off +Notes Audit Transfusion Notes Favourites

- Ability to 'sign off' reports
- Ability to add notes
- Access the audit for the selected UR number
- Access the transfusion history, available products and antibody register
- Access to Notes
- Access to Favourites

My Tasks Enquiry User Settings

- Access Enquiry screen to search another patient
- Access My Tasks and User Settings

Displays the next report to be signed off

Next Unsigned Report

Results across linked UR numbers, DO NOT display on Episodes view

Log out and external links

Episode Time and date is Collection time and date

- UE
- 09:10 03-May-2017
- CK, LD, UA
- 09:05 03-May-2017
- CMP
- 09:30 03-May-2017
- ULB
- 06:03 03-May-2017
- AMY, LIPASE
- TROP
- 06:00 03-May-2017
- CK, LD, UA
- CMP, ULB
- 08:30 02-May-2017
- ITD
- TR
- ALB, AMY, AST, BILC, CA, CHOL, CK, CRP, LIPASE, UA

Icon status displayed

Arrow shows episode displayed

Orange Text shows results displayed on report

auscare KILLINGTON

Investigations Ordering Sign Off +Notes Audit Transfusion Notes Favourites

Dummy, Patient Alfred PM981109999

HCF: CLIN UNIT: DOB: 12-Nov-1936 SEX: Male

WARD: UNKNOWN DOCTOR: CONSULTANT

Peter Mac Pathology 305 Grattan St, Melbourne (03) 8559 5000

Accredited for Compliance with NPAAC Standards and ISO 15189

Report Delivery Address: PM 2C - 2D Specialist Clinics

Lab No: 80016272 Collected: 03-May-17 09:10

Received: 03-May-17 09:11 Validated: 02-Jun-17

Ext Ref: Trial Ref: Specimen: BLOOD

Requester: Unknown Doctor

See below for copies

View Request form by selecting the lab number

UR: PM981109999 DOB: 12-Nov-1936 Age: 80 years SEX: M Address: 99 MOLLUSC AVENUE CARLTON Address: 96561234 Phone: 96561234 Location: Peter MacCallum Cancer Centre

Most recent

Oldest results

Units and Reference Ranges

Laboratory:	PM 03-May-17	Units	Ref. Interval				
Time:	06:00	06:03	09:00	09:05	09:10		
Episode No.:	80016266	80016266	80016266	80016272	80016272		
Spec. Type:	BLOOD	BLOOD	BLOOD	BLOOD	BLOOD		
General Chemistry							
Sodium	137		140			mmol/L	(135 - 145)
Potassium	4.2		3.9			mmol/L	(3.5 - 5.2)
Chloride	106		99			mmol/L	(95 - 110)
Bicarbonate	23		27			mmol/L	(22 - 32)
Urea	7.0		10.8 H			mmol/L	(3.0 - 9.2)
Creatinine	97		224 H			umol/L	(60 - 110)
ALP	63 L		23 L				See Below (> 90)
hs Troponin I		10				ng/L	(< 26)
Creatine Kinase	78					U/L	(0 - 200)
Calcium							
Calcium	2.19		2.21		2.90 H	mmol/L	(2.10 - 2.80)
Calcium Corrected	2.09 L		2.19		3.00 H	mmol/L	(2.10 - 2.60)
Magnesium	0.79		1.16 H		0.82	mmol/L	(0.70 - 1.10)
Phosphata	1.77 H		1.55 H		1.16	mmol/L	(0.75 - 1.50)
Total Protein							
Albumin	45		41		35	g/L	(35 - 50)
Globulin	21 L		19 L			g/L	(25 - 40)
Bilirubin	< 2		20			umol/L	(< 21)
ALT	21		291 H			U/L	(5 - 40)
SGT	30		128 H			U/L	(< 60)
ALP	70		67			U/L	(90 - 110)
LD	241				348 H	U/L	(120 - 250)
Lipase							
		24				U/L	(< 60)

View Graph by selecting the test name and selecting Graph button **Graph**

Print result

Scroll bar

ICONS: A coloured page or an hourglass will display if there are unviewed/unsigned off/pending pathology results in the list of results using Reports, Requests or Episodes view types.



- If Status has exceeded the configured '**overdue**' time the icons appear as **flashing** and will remain **flashing** until they are signed off

This function enables Clinicians to 'sign off' results have been viewed. This action changes the icon colour.

This function enables Notes to be added to the report being viewed. These notes are audited.

View audit

View Notes

This function enables Clinicians to view the next result that has not been signed off.

Select Favourites for the drop down menu options to display.

Favourites

- + Favourite ... Select List
- Favourite ... Personal-Haem Clinic

Patients may be added or removed from Favourite Lists which are created via User Settings

The Transfusion button displays when the patient has crossmatch or products assigned.

Product Availability		Antibody Register	
Transfusion History			
Unit No	Product	Group	Status
8777777	RCC SAGM IRR	O NEG	Compatible

Select Transfusion to display the Transfusion History screen and access buttons for Product Availability screen and Antibody Register screen.

Back Use the Back button to return to the previous screen

T/F = To Follow
Interim = Pending